Respiratory protection Options & Current Standards Guidance

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SEA OR&CSSD Scientific Affairs & Education Manager

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Procedure Mask

- Disposable, simple mask used to cover nose and mouth

Surgical Mask

- Disposable, mask used to cover nose and mouth
- Provides a fluid barrier

Filtering Facepiece Respirators (FFRs), such as N95

- Disposable, respirator used to cover nose and mouth
- Fluid resistant (Surgical N 95 Respirator)

Elastomeric Respirator

- Reusable respirator
- **Half & full** facepiece offerings
- Multi-piece: facepiece and **disposable filters**

Powered Air Purifying Respirator (PAPR)

- Reusable respirator
- Multiple face/head cover offerings
- Multi-piece: head cover, tube, battery pack and belt

https://www.3m.com/3M/en_US/worker-health-safety-us/covid19/

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Masks are not designed to seal to the face.

Surgical Masks

help prevent particles expelled by wearers from entering sterile field.

Surgical Masks

help prevent high-velocity sprays from reaching the healthcare professional's mouth and nose.

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
American Society of Testing and Materials(ASTM) F2100-11

1. BFE (Bacterial Filtration Efficiency) -ASTM F2101
2. PFE (Particulate Filtration Efficiency) -ASTM F2299
3. Fluid Resistance – ASTM F1862
4. Delta P (Pressure Differential)
5. Flame Spread – 16 CFR Part 1610

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The US Center for Disease Control (CDC) cites the N95 respirator standard as part of the advised protective equipment in their [Covid-19 FAQ](#) and their [SARS guidance](#) (SARS being a similar type of Corona virus). Which suggests that an **N95 or better respirator is acceptable.**
Crisis : equivalence to N95 from other countries

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


Respiratory Protection.

How N95s work

- FFRs contain many layers of crisscrossing, non-woven, fibrous filter material that capture particles
- As the airstream flows through the layers of filter media, particles are trapped/captured, by several mechanisms:
 - Very large particles in the airstream will settle out **due to gravity**
 - Other particles **may impact a fiber and be captured**
 - Very small particles are captured by **diffusion**

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**NIOSH Certified
or Approved*
TC number**

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NIOSH Certified
N95 Respirators
TC –Number

*Surgical N95
Respirators*

Surgical Masks

FDA Cleared

* <https://wwwn.cdc.gov/niosh-cel/>

https://www.cdc.gov/niosh/npptl/topics/respirators/disp_part/

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Limited Reuse

Non NIOSH approved
beyond shelf life

Decontamination
& Reuse (except AGP)

<https://www.cdc.gov/coronavirus/2019-ncov/hcp/respirators-strategy/index.html> Accessed April 12, 2020

N95 Respirator.

Decontamination (last resort)

3M is collaborating with several sterilization companies and institutions that are investigating ways for hospitals to safely decontaminate 3M's N95 FFRs in line with the CDC guidance on Decontamination and Reuse of Filtering Facepiece Respirators.

<https://www.cdc.gov/coronavirus/2019-ncov/hcp/ppe-strategy/decontamination-reuse-respirators.html>

Four key aspects of successful decontamination reprocessing:

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Per the CDC guidelines, a number of sterilization companies are assessing **decontamination processes for N95 filtering facepiece respirators (FFRs).**

The U.S. Food and Drug Administration (FDA) is evaluating and granting Emergency Use Authorizations (**EUAs**) for sterilization companies assessing decontamination processes for N95 filtering facepiece respirators (FFRs) during the COVID-19 outbreak.

Issued EUAs for Personal Protective Equipment with regards to COVID-19 will be available on the **FDA website: Personal Protective Equipment EUAs**

<https://multimedia.3m.com/mws/media/18248690/decontamination-methods-for-3m-n95-respirators-technical-bulletin.pdf>
<https://www.cdc.gov/coronavirus/2019-ncov/hcp/ppe-strategy/decontamination-reuse-respirators.html>

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Technical Bulletin

May, 2020
Revision 6

Decontamination Methods for 3M Filtering Facepiece Respirators such as N95 Respirators

Background

NOTE: Please revisit this document often for frequent updates.

During this public health emergency of the COVID-19 pandemic outbreak, many healthcare institutions are experiencing shortages of filtering facepiece respirators such as N95 respirators. The U.S. Center for Disease Control and Prevention (CDC) has issued [Strategies for Optimizing the Supply of N95 Respirators](#). In this document the CDC recommends conventional capacity strategies, contingency capacity strategies (during expected shortages) and crisis strategies (during known shortages). Contingency and crisis strategies include use of N95s past their shelf life, extended use of N95s, use of other types of respirators, use of respirators from other countries, and re-use of respirators, ahead of decontamination of respirators.

The CDC discusses reuse and extended use of N95s as a Crisis strategy at [Recommended Guidance for Extended Use and Limited Reuse of N95 Filtering Facepiece Respirators in Healthcare Settings](#) and has published guidelines on [Decontamination and Reuse of Filtering Facepiece Respirators](#). CDC says research indicates the virus survives for up to 72 hours on a variety of surfaces. Therefore, CDC is recommending a wait and reuse approach before consideration of other decontamination approaches.

<https://multimedia.3m.com/mws/media/1824869O/decontamination-methods-for-3m-filtering-facepiece-respirators-technical-bulletin.pdf>

Table 1: Effect of decontamination methods on certain 3M N95 Filtering Facepiece Particulate Respirators (Continued) (Sheet 2 of 2)

Decontamination Method	3M N95 Models Evaluated ^a	Cycle	Number of Reprocess Cycles Tested	Filtration Efficiency ^b	Fit Related Evaluation	U.S. FDA EUA Issued
VHP- Battelle	1860, 8210, 1804	Under evaluation	3 cycles: tested 20 cycles: Under evaluation	3 cycles: Pass 20 cycles: Under evaluation	3 cycles: Pass 20 cycles: Under evaluation	EUA
						Facility Instructions
						HC Personnel Instructions
VHP – Steris - Victory™, 1000 ED, ARD, and M100 Biodecontamination Unit	1860, 8210, 1804, 1870+	STERIS Atmospheric VHP Process	20	Pass	Pass	No
Ultraviolet Light Environmental Decontamination Systems						
UV-C (254 nm)	1860, 8210, 1804	Refer to CDC guidance ^d or UV OEM	Maximum 100 J/cm ² cumulative lifetime exposure	Pass	Pass	No
Xenex Lightstrike™ System	1860, 8210, 1804	Pulsed xenon, 200 – 280 nm for 5 minutes	10	Pass	Pass	No

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Method	Treatment level	FFR filtration performance	FFR fit performance	Other observations	References
Vaporous hydrogen peroxide (VHP)	Battelle report: Bioquell Clarus C HPV generator: The HPV cycle included a 10 min conditioning phase, 20 min gassing phase at 2 g/min, 150 min dwell phase at 0.5 g/min, and 300 min of aeration.	Passed	FFR fit was shown to be unaffected for up to 20 VHP treatments cycles using a head form	Degradation of straps after 30 cycles (Battelle report)	3, 4
Ultraviolet germicidal irradiation (UVGI)	0.5–950 J/cm ²	Passed	90–100% passing rate after 3 cycles depending on model		2, 3, 7, 8, 9, 10
Moist heat incubation	15–30 min (60°C, 80% RH)		H1N1 influenza A/PR/8/34	99.99%	14

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